


# Wurapa Small Joint Compressor and Distractor

Designed by Raymond K. Wurapa, MD

*Designed to allow one-handed manipulation and deployment once fixation pins are placed*



Now available with two hole sizes on each instrument!



Designed to simplify several small joint procedures:

- ▶ Preparation of small bone non-unions before bone grafting and fixation
- ▶ Preparation of small joints for arthrodesis (e.g. partial wrist fusion)
- ▶ Useful for intercarpal stabilization while performing ligament reconstructions (e.g. scapholunate ligament repair/reconstruction)
- ▶ Distract and better evaluate small joints before determining final management.

Pins should be cut short above the pin guides to allow full access to the operative site.

PRODUCT NO.'S:	
Double .045" (1,1 mm) & .062" (1,6 mm) Holes	Single .045" (1,1 mm) Hole
1751 [Compressor] Compresses From: 28 mm Overall Length: 4.625" (11,7 cm)	1753 [Compressor] Compresses From: 28 mm Overall Length: 4.5" (11,4 cm)
1752* [Distractor] Distracts to: 46 mm Overall Length: 4.625" (11,7 cm)	1754 [Distractor] Distracts to: 46 mm Overall Length: 4.5" (11,4 cm)

MADE EXCLUSIVELY FOR INNOVIMED IN GERMANY

\* PROUDLY MADE IN THE USA

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to Launch Our Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**