


Wurapa Swivel Blade Forearm Retractor

Designed by Raymond Wurapa, MD

Designed for forearm and wrist fracture exposure, the blades swivel for less stress on soft tissue

Swivel-blade technology helps to allow parallel deployment of retractor blades to maximize wound exposure and minimize edge loading on surrounding soft tissues. Parallel deployment of the retractor blades also helps prevent rotation and migration of the retractor during a procedure.



PRODUCT NO'S:

1646-00 [Set]
Includes Retractor and Two Swivel Blades

Also available individually:

1646-01 [Retractor]
Overall Length: 5.125" (13 cm)
Opens to: 2.5" (6,4 cm)

1646-02 [Swivel Blade]
One blade with this product number, two included in set
Width: .9375" (24 mm)
Depth: .75" (19 mm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2017 Innomed, Inc.