



# Extra Large Hip Retractors

Designed by Wayne M. Goldstein, MD

*For hip surgery with large patients, and when extra large instruments are desired for increased leverage and depth*


## Extra Leverage Femoral Neck Elevators

PRODUCT NO'S:	
7650 [Standard]	Overall Length: 18.25" (46,4 cm) Handle Length: 9.25" (23,5 cm) Blade Width: 38mm
7650-02 [Short Handle]	Overall Length: 15.25" (38,8 cm) Handle Length: 6.25" (15,9 cm) Blade Width: 38 mm



## Infero-posterior Acetabular Capsule Retractors

PRODUCT NO'S:	
7620-01 [Right]	Overall Length: 12" (30,5 cm) Handle-to-Bend Length: 6" (15,2 cm)
7620-02 [Left]	Overall Length: 12" (30,5 cm) Handle-to-Bend Length: 6" (15,2 cm)



## Extra Leverage Proximal Femoral Elevator

PRODUCT NO:	
7640	Overall Length: 17.5" (44,5 cm) Handle Length: 13" (33 cm) Blade Width at Widest: 63 mm



## Large Cobra Retractors

PRODUCT NO'S:	
7630-01 [Standard]	Overall Length: 17.5" (44,5 cm) Handle Length: 14" (35,6 cm)
7630-02 [Wide]	Overall Length: 17.5" (44,5 cm) Handle Length: 14" (35,6 cm)



## Modular Weights

*Weights can be used to help hold the retractors in place.*

PRODUCT NO'S:	
3430-01	1.5 lbs. (.68 kg)
3430-02	2.0 lbs. (.91 kg)
3430-03	2.5 lbs. (1.13 kg) with attaching hook



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**