




Extra Long Rongeur

Helpful in minimally invasive total hip surgery by keeping hands out of the field of view



Rongeur Jaws at Actual Size



PRODUCT NO:
1771-01 Jaw Bite: 5 x 16 mm Overall Length: 14" (35,6 cm)
1771-02 Jaw Bite: 6 x 16 mm Overall Length: 14" (35,6 cm)
1771-03 Jaw Bite: 12 x 16 mm Overall Length: 14" (35,6 cm)



Mazzara Pistol Grip Extra Long Rongeur

Designed by James T. Mazzara, MD

Pistol Grip handle lessens hand fatigue and slippage, and allows for better visualization

PRODUCT NO:
1768-02 Jaw Bite: 7 x 16 mm Overall Length: 12.5" (31,7 cm) Shaft-to-End Length: 6" (15,2 cm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.