

# Zelicof Winged Retractor

Designed by Steven B. Zelicof, MD Ph.D


*A non-modular, fixed point retractor for pelvic visualization in THA and pelvic surgery, with the distal post and wing providing bony stability and retraction*



Thin Shaft  
#6117-5

Thick Shaft  
#6117-10

PRODUCT NO:	
6117-5	[Thin Shaft]
Shaft Diameter: 6 mm	
Overall Length: 9" (22.9 cm)	
Pad Dimensions: 31,4 x 19,6 mm	
Prong Length: .85" (21,6 mm)	
Prong Width (before taper): 3 mm	
6117-10	[Thick Shaft]
Shaft Diameter: 10 mm	
Overall Length: 9" (22,9 cm)	
Pad Dimensions: 31,4 x 19,6 mm	
Prong Length: .85" (21,6 mm)	
Prong Width (before taper): 3 mm	



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

1.800.548.2362